




PHOTOBIOREACTOR

A photobioreactor uses light to cultivate photosynthetically active microorganisms such as algae, cyanobacteria or phototrophic bacteria. Of particular note is the advanced lighting supply and in-situ measurement technology for monitoring light irradiation in the bioreactor. For this purpose, LED light bars and an in-situ multi-channel sensor system are used.

Components overview



- ① Gas quality analyzer
- ② Mixer unit
- ③ Photobioreactor
- ④ Lighting system
- ⑤ Heating / cooling system
- ⑥ Sparkling system
- ⑦ Gas amount meter
- ⑧ SENSOcontrol data acquisition and control system

The photobioreactor provides controlled conditions for:


- ✓ Light supply
- ✓ Temperature
- ✓ Nutrient supply
- ✓ Gas management

Features and options

- ✓ Up to 10 bioreactors on one gas quality measurement system
- ✓ Reactor with different working volumes
- ✓ Continuous operation with feeding
- ✓ Electrical heating / cooling
- ✓ SENSOcontrol with touch panel
- ✓ Sparkling system with MFC (mass flow controller)
- ✓ Dosage pumps for pH-control
- ✓ Advanced lighting system
- ✓ Measurement equipment: Gas quantity, pH, redox, temperature, gas quality and more




PHOTOBIOREACTOR




SENSOcontrol

The universal automation system SENSOcontrol is used for data acquisition applications and control of units (pumps, mixer) of the test plant. A touch panel for intuitive operation is integrated for user-friendly parameter set-up and real time visualization of measured data.



In-situ light intensity sensor

A light intensity sensor with 14 channels is available. This allows semi-quantitative measurement of in-situ light intensity. The sensor measures wavelengths between 350 and 900 nm.



Sparkling systems

The additional sparkling system with a mass flow controller allows the defined sparkling of gases. Gas mixing systems are also available.



LED lighting system

- ✓ **Red wavelength:** 670 nm, each with 15 LEDs; total power per module up to 1.5 W
- ✓ **IR wavelength:** 740 nm, each with 15 LEDs with integrated visible control LEDs; total power per module up to 1.5 W
- ✓ **Green wavelength:** 570 nm, each with 15 LEDs; total output per module up to 1,200 lm/2 W
- ✓ **Blue wavelength:** 480 nm, each with 15 LEDs; total output per module up to 270 lm/3.6 W
- ✓ **White color:** 2,700 K with 15 LEDs each; total output per module up to 825 lm
- ✓ Dimensions of each lighting module 300 x 15 mm – others on request

Information on Photobioreactor



The equipment features will be defined with the quotation.

© 2025 UIT | All rights reserved | V1